

MICROLUTE® CP - SAX

PROTOCOL


Example Method for the Extraction of Weak Acids

The Microlute® CP strong anion exchange (SAX) uses a quaternary ammonium chemistry immobilised on a polymeric base with a pK_a of >18, making it ideal for the capture of weak acidic analytes through anion exchange. As with all Microlute® CP products, the polymeric base offers a secondary retention mechanism for neutral compounds.

	Solvent	10 mg Volume	2 mg Volume
1 Condition	MeOH	500 µL	200 µL
2 Equilibrate	0.1% ammonium hydroxide in H ₂ O	500 µL	200 µL
3 Load	Sample diluted with 0.1% ammonium hydroxide in H ₂ O	500 µL*	100 µL*
4 Wash 1	0.1% ammonium hydroxide in H ₂ O	500 µL	200 µL
5 Wash 2	0.1% ammonium hydroxide in MeOH	500 µL	200 µL
6 Elute	2% formic acid in MeOH	500 µL	2x 50 µL
7 Analyse	Dilute eluent, directly inject or evaporate eluent and reconstitute in a more suitable solution for analysis.		

This method is an ideal starting point for several applications and for samples containing a wide range of components. Method development may be required to obtain optimal recovery and reproducibility.

*The wells can hold up to 1 mL for the 2 mg product and 1.5 mL for the 10 mg product. However, the actual volume of the sample that can be loaded is determined by its concentration and how much it has been diluted.



porvair
sciences

Technical Support & Sales

Tel: +44 1978 661144

Technical: technical@porvairsciences.com

Sales: int.sales@porvairsciences.com

Learn more

www.microplates.com

Porvair and Microlute® are registered trademarks of Porvair plc.

© Copyright 2025. Porvair Sciences Ltd. All rights reserved.

While every effort has been made to ensure the accuracy of this document, due to continuous product development, the data contained is subject to constant revision and Porvair Sciences reserves the right to change, alter or modify its contents.