



REV.	ECO NO.	DESCRIPTION
A	000002	RELEASE TO PRODUCTION.



- NOTES:  
 1. MAX VOLUME PER WELL: 815  $\mu$ l  
 2. SURFACE AREA OF FILTER: 19.35  $\text{mm}^2$

CATALOG NO.	FILTER
F20002	UHMW PE 25 $\mu$ m
F20003	GLASS FIBER 1.2 $\mu$ m
F20004	PVDF 0.45 $\mu$ m
F20005	PP 0.45 $\mu$ m
F20007	GLASS FIBER 1.0 $\mu$ m
F20008	GLASS FIBER 0.7 $\mu$ m
F20056	PC 2.0 $\mu$ m
F20068	PE FRIT + GLASS FIBER 0.7 $\mu$ m



**INNOVATIVE**  
 CUSTOM MICROPLATE  
 MICROPLATES

1998 WESTOVER RD.  
 CHICOPEE, MA 01022  
 TEL. (413) 593-5300  
 FAX (413) 593-5900

© 2001  
 CONFIDENTIAL  
 TRADE SECRETS  
 PATENTS ISSUED  
 AND PENDING  
 DO NOT COPY  
 OR DISCLOSE  
 WITHOUT WRITTEN  
 PERMISSION OF  
 INNOVATIVE MICROPLATE

TOLERANCES UNLESS  
 OTHERWISE SPECIFIED  
 XX  $\pm$  .03  
 XXX  $\pm$  .015  
 XXXX  $\pm$  .0005  
 FRACTION:  $\pm$  1/64  
 ANGULAR:  $\pm$  0°15'  
 RADIUS 0.015" UNLESS  
 OTHERWISE SPECIFIED  
 BREAK ALL SHARP EDGES

**MATERIAL:** RESPPCP  
 NATURAL POLYPROPYLENE

**CONCENTRATE:**

**DATE:** 08/10/01  
**DRAWN BY:** YELENA SUKONIK

**PRT NAME:**  
 96 WELL, FILTER PLATE, 800 $\mu$ l, LONG DRIP

**CATALOG NO.:** SEE ABOVE

<b>APPR BY:</b>	<b>PRT SIZE:</b> A	<b>SCALE:</b> NTS	<b>SHEET</b> 1 OF 1	<b>REV.</b> A
-----------------	-----------------------	----------------------	---------------------	------------------