




REV.	ECO NO.	DESCRIPTION
A	000022	RELEASE TO PRODUCTION.



NOTE:

MAX VOLUME OF ROW 32 ml

<b>INNOVATIVE</b> CUSTOM MICROPLATE MICROPLATES 1998 WESTOVER RD. CHICOPEE, MA 01022 TEL. (413) 593-5300 FAX (413) 593-5900	© 2001 CONFIDENTIAL TRADE SECRETS PATENTS ISSUED AND PENDING DO NOT COPY OR DISCLOSE WITHOUT WRITTEN PERMISSION OF INNOVATIVE MICROPLATE	TOLERANCES UNLESS OTHERWISE SPECIFIED XX ± .03 XXX ± .015 XXXX ± .0005 FRACTION: ± 1/64 ANGULAR: ± 0°15' RADIUS 0.015" UNLESS OTHERWISE SPECIFIED	<b>MATERIAL:</b> RESPPCP NATURAL POLYPROPYLENE <b>CONCENTRATE:</b>	<b>PRT NAME:</b> 8 ROW, PARTITIONED RESERVOIR, PYRAMID BOTTOM, NATURAL PP	<b>SCALE:</b> NTS	<b>SHEET</b> 1 OF 1	<b>REV.</b> A
	<b>DATE:</b> 01/22/02 <b>DRAWN BY:</b> YELENA SUKONIK	<b>PRT NO.:</b> S30020	<b>APPR BY:</b> A	<b>PRT SIZE:</b> A	<b>MAX VOLUME OF ROW</b> 32 ml	<b>REV.</b> A	