



REV.	ECO NO.	DESCRIPTION
A	000316	RELEASE TO PRODUCTION.



NOTE:
1.) VOLUME RESERVOIR 170ml

FOR REFERENCE ONLY

Seahorse Bioscience
300 Griffith Road
Chicopee, MA 01022
Tel: (413) 593-5300
Fax: (413) 593-5900

© 2001 CONFIDENTIAL TRADE SECRETS PATENTS ISSUED AND PENDING DO NOT COPY OR DISCLOSE WITHOUT WRITTEN PERMISSION OF INNOVATIVE MICROPLATE	TOLERANCES UNLESS OTHERWISE SPECIFIED XX ± .03 XXX ± .015 XXXX ± .0005 FRACTION: ± 1/64 ANGULAR: ±0°15' RADIUS 0.015" UNLESS OTHERWISE SPECIFIED BREAK ALL SHARP EDGES	MATERIAL: RESPPCP NATURAL POLYPROPYLENE		PRT NAME: 384 WELL RESERVOIR WITH PYRAMID BOTTOM				
		CONCENTRATE:		CATALOG NO.: S30042				
		DATE: 5/19/05	DRAWN BY: KN	APPR BY:	PRT SIZE: A	SCALE: NTS	SHEET 1 of 1	REV. A