


96-well Microplate (PP) Deep Well (Square)



L	127.60 ± 0.10	H350	15.30 ± 0.10
W	85.50 ± 0.10	H1	24.50 ± 0.10
H	44.40 ± 0.05	H2	39.90 ± 0.05
Hs	4.00 ± 0.10	P	9.00 ± 0.05
X1	14.00 ± 0.05	D	8.25 ± 0.05
Xn	113.50 ± 0.05	C	12
Y1	11.00 ± 0.05	R	8
Yn	74.50 ± 0.05		

Volume 350µl = 750µl, 1ml = 1,400µl, 2ml = 2,500µl
 Well form V Bottom
 Material PP