


## Component Schematic

**Krystal 96 Microplate**

**P/N 214nnn and 215nnn**



L	127.60 ± 0.25	X1	14.40 ± 0.10	C	12
W	85.60 ± 0.25	Xn	113.40 ± 0.10	R	8
H	14.70 ± 0.10	Y1	11.50 ± 0.10	P	9.00 ± 0.10
Hs	6.00 ± 0.10	Yn	74.50 ± 0.10	D	6.30 ± 0.10
Hw	11.15 ± 0.10	Hc	2.75 ± 0.10	after assembly	

Dimensions in mm

Well volume 350µl, after assembly

Well form Flat bottom

Material GPPS